



Cleanliness



CC G4

User-friendly
cleanliness cabinet for
accurate particle extraction

Specialized - User friendly - Ergonomic



Features:

- Timesaving

- By blank tests
- By extraction

- Ergonomic



- Height adjustable for operator
- Control rinsing flow by foot pedal

- Cascade filter mount


- Cartridge for both D47 and D90 filters

Functions:

- Touchscreen
- Flow control in bar or lpm
- Preset amount of rinsing fluid with automatic stop
- Reminder for service intervals
- Settings for either alkaline chem or degreaser solution



Complete provider of laboratory consumables



References

VOLVO

Construction Equipment



Donaldson®
Filtration Solutions

Hydra
Specma



BOSCH

rexroth

A Bosch Company



Cleanliness Cabinet CC-G4

Special designed cabinet for particle extraction of industrial components according to ISO16232, ISO4406/4407, VDA19.1 and custom standards.

Developed with 50 years of expertise in metrology and 20+ years of focus in technical Cleanliness.

Why choose NYLI

- Design, development, assembly and quality check all made in our facilities in Sweden
- Installation, service, training and support by specialized engineers trained by NYLI
- NYLI have an accredited laboratory for particle testing working with customers cleanliness specifications daily in our own systems.
- CE and RoHS compatible

